

#### SCOPE OF ACCREDITATION TO ISO/IEC 17025:2017 & ANSI/NCSL Z540-1-1994

## METTLER TOLEDO, LLC 2915 Argentia Road, Unit 6

Mississauga, Ontario, CANADA L5N 8G6 Tamara Rogers-McCrossan Phone: +1-416-708-7781

#### **CALIBRATION**

Valid To: January 31, 2023 Certificate Number: 1788.02

In recognition of the successful completion of the A2LA evaluation process, accreditation is granted to this laboratory to perform the following calibrations<sup>1</sup>:

#### I. Fluid Quantities

Parameter/Equipment	Range	CMC <sup>2</sup> (±)	Comments
Piston Operated Volumetric Apparatuses (POVA): Piston Pipettes, Piston Burettes, Dilutors & Dispensers <sup>3</sup>	Up to 2 μL (>2 to 10) μL (>10 to 20) μL (>20 to 100) μL (>100 to 200) μL (>200 to 1000) μL (>200 to 1000) μL (>1000 to 2000) μL (>2000 to 5000) μL (>5000 to 10 000) μL (>10 000 to 20 000) μL (>10 000 to 20 000) μL (>10 000 to 20 000) μL	0.016 μL 0.022 μL 0.044 μL 0.27 μL 0.38 μL 2.2 μL 3.5 μL 7.1 μL 14 μL 23 μL	Gravimetric method per ISO-8655-6

<sup>&</sup>lt;sup>1</sup> This laboratory offers commercial calibration service.

Page 1 of 1

<sup>&</sup>lt;sup>2</sup> Calibration and Measurement Capability Uncertainty (CMC) is the smallest uncertainty of measurement that a laboratory can achieve within its scope of accreditation when performing more or less routine calibrations of nearly ideal measurement standards or nearly ideal measuring equipment. CMCs represent expanded uncertainties expressed at approximately the 95 % level of confidence, usually using a coverage factor of k = 2. The actual measurement uncertainty of a specific calibration performed by the laboratory may be greater than the CMC due to the behavior of the customer's device and to influences from the circumstances of the specific calibration.

<sup>&</sup>lt;sup>3</sup> This laboratory performs field calibration activities for these parameters.



# **Accredited Laboratory**

A2LA has accredited

## METTLER TOLEDO, LLC

Mississauga, Ontario, CANADA

for technical competence in the field of

### Calibration

This laboratory is accredited in accordance with the recognized International Standard ISO/IEC 17025:2017

General requirements for the competence of testing and calibration laboratories. This laboratory also meets ANSI/NCSL Z540-1-1994 and R205 – Specific Requirements: Calibration Laboratory Accreditation Program. This accreditation demonstrates technical competence for a defined scope and the operation of a laboratory quality management system (refer to joint ISO-ILAC-IAF Communiqué dated April 2017).

SEAL 1978 V JOHN A2LA

Presented this 25th day of January 2021.

Vice President, Accreditation Services For the Accreditation Council

Certificate Number 1788.02

Valid to January 31, 2023